

Virus Transport Medium (VTM)

Country of Origin :		
Name of Manufacturer :		
Brand Name :		
Item name	Virus Transport Medium (VTM)	Bidder's Offer
Purpose	Swab collection (COVID 19)	
Tube Requirements	Suitably prepared sterile media for use in collecting throat and nasal swabs from human. Should contain virus inactivator, should stabilize virus RNA. Should be contained in airtight plastic tubes with cap. There must be sticker for labeling	
Media volume	1ml to 3ml	
Swab collection sticks	<ul style="list-style-type: none"> - Along with the tubes there should be two swab sticks (one for oropharyngeal swab and another for nasopharyngeal swab). - There should be provision of break lines to allow to fit into the container. - Both the sticks should have fiber acrylic swab. The stick for nasopharyngeal swab should be flexible enough for ease the collection of swab sample from nasopharynx. - The swabs sticks should be contained in blisters 	
Others	The item should be atleast CE approved	